



# Specification

#### General

The Powerex medical vacuum system is designed to create a suction system to remove unwanted fluids or gases from hospital/ laboratory working areas. The medical vacuum system package is compliant with the NFPA 99 requirements for Risk Category 1 systems. Each system is completely tested before shipment and

#### includes:

- Multiple vacuum pumps and associated equipment.
- · AMSE air receiver.
- Medical control panel.

Each pump is factory piped to a common intake manifold. Vibration isolation pads are included with the system.

### **Claw Vacuum Pump**

Each pump shall be a rotary claw type vacuum pump, and shall be direct-driven through a shaft coupling by a C-face, TEFC electric motor.

- Each vacuum pump shall be dry-running, featuring two claw-type, non-contacting rotors and shall not require any sealing fluid in the pumping chamber, assuring virtually maintenance-free operation.
- Each vacuum pump shall include an internal relief valve, and a built-in, anti-suck-back valve mounted at the pump inlet.
- Each pump within the system shall include a check valve, inlet and discharge flex connectors, a 5 micron inlet filter and a pump isolation valve.

# Motor

The motor is continuous duty, C-face, TEFC, suitable for 208-230, or 460V, 3 phase, 60 hertz electrical operation.

#### Air Receiver

The system shall include an ASME rated air receiver. The tank shall be equipped with a vacuum gauge, a sight gauge, by-pass valves, and a manual drain.

### **Premium NFPA Control Panel**

The control system provides automatic lead/lag sequencing and automatic alternation of all pumps in order to equalize the amount of usage among the available vacuum pumps. The Premium NFPA Control panel shall include a gateway server card and all features listed below:

- PLC controller and a color touch screen panel which displays the operating status of the unit.
- Building automation communication gateway with BacNet® protocol and Web server features. Web servers features include email notifications in case the system is in alarm or has achieved one its maintenance intervals and requires service.
- Ethernet port for connection to BacNet® server or direct connection to facility Ethernet for viewing of system operations and status via device IP-address.
- UL508A listed control panel in a NEMA 12 enclosure. The panel door will include: the HMI touch screen, an audible and visual

alarms with an acknowledge button, and an HOA switch for each pump.

- Magnetic starters.
- Vacuum transducer for process control.
- Single point power connection.
- Redundant 120Vac control transformers with fused primary and secondary protection.
- System overload trip, high temperature conditions or maintenance intervals for the pump will result in visual and audible alarms.

### **Optional VFD Control Panel**

Variable Frequency Drive (VFD) control improves efficiency over a conventional "on/off" demand based system by more closely matching the pump speed to the changing load requirements. All VFD systems come standard with a Premium NFPA Control panel. The VFD and Premium NFPA Control Panel shall include a gateway server card and all features listed below:


- PLC controller and a color touch screen panel which displays the operating status of the unit.
- Building automation communication gateway with BacNet® protocol and Web server features. Web server s features include email notifications in case the system is in alarm or has achieved one its maintenance intervals and requires service.
- Ethernet port for connection to BacNet® server or direct connection to facility Ethernet for viewing of system operations and status via device IP-address.
- UL508A listed control panel in a NEMA 12 enclosure. The panel door will include: the HMI touch screen, an audible and visual alarms with an acknowledge button, VFD start/stop switch and an HOA switch for each pump.
- · Magnetic starters.
- Vacuum transducer for process control.
- Single point power connection.
- Redundant 120Vac control transformers with fused primary and secondary protection.
- System overload trip, high temperature conditions or maintenance intervals for the pump will result in visual and audible alarms.


### **Available Options**


- Oxygen assured pumps prepared for use in WAGD systems.
- Internal tank lining for corrosion resistance.
- Variable Frequency Drive (VFD) Control Panel with Premium NFPA controls.



Dimensions										
Model	Dim. A	Dim. B Dim. C		Inlet	Outlet (x4)					
CVTD0203V	35"	56"	82"	1 ½"	1"					
CVTD0303V	35"	56"	82"	1 ½"	1"					
CVTD0504AV	35"	56"	85"	2"	1"					
CVTD0504BV	35"	56"	85"	2"	1"					







Claw Vacuum Tankmount Duplex with Premium Controls											
Model HF	LID S	SCFM @ 19"	NFPA System	Tank Size	BTU/Hr (2)	DTI / Ur (2) dB(A)	System F.L.A.			System	
	ПР	Hg each	Capacity <sup>(1)</sup>	(gal)		Level (3)	208V	230V	460V	Weight (lbs)	
CVTD0203V	2 (2)	16	16	80 V	4,335	70	14.4	13.2	7.6	1,300	
CVTD0303V	3 (2)	21	21	80 V	6,503	70	20.8	18.8	10.4	1,520	
CVTD0504AV	4 (2)	29	29	120 V	8,670	79	22	20.4	11.2	1,650	
CVTD0504BV	5 (2)	38	38	120 V	10,838	79	26.8	24.8	13.4	1,720	

## Notes:

- 1 System Capacity is shown with one or more pumps in reserve per NFPA 99.
- 2 BTU/Hr levels are shown with reserve pump(s) on standby.
- 3 dB(A) levels are shown with one pump in reserve per NFPA 99.
- 4 3 Year Limited Warranty.