

Specification

General

The Powerex medical vacuum system is designed to create a suction system to remove unwanted fluids or gases from hospital/laboratory working areas. The medical vacuum system package is compliant with the NFPA 99 requirements for Risk Category 1 systems. Each system is completely tested before shipment and includes:

- Multiple vacuum pumps and associated equipment.
- AMSE air receiver.
- Medical control panel.

Each pump is factory piped to a common intake manifold. Vibration isolation pads are included with the system.

Oilless Dry Rotary Vane Vacuum Pump

Each pump shall be a dry rotary vane type vacuum pump, and shall be direct-driven through a shaft coupling by a C-face, TEFC electric motor.

- Each vacuum pump shall be dry-running, featuring self-lubricating carbon/graphite vanes and shall not require any sealing fluid in the pumping chamber.
- Each vacuum pump shall include an internal relief valve, a check valve, inlet and discharge flex connectors, a 5 micron inlet filter and a pump isolation valve.

Motor

The motor is continuous duty, C-face, TEFC, suitable for 208-230, or 460V, 3 phase, 60 hertz electrical operation.

Air Receiver

The system shall include an ASME rated air receiver. The tank shall be equipped with a vacuum gauge, a sight gauge, by-pass valves, and a manual drain.

Premium NFPA Control Panel

The control system provides automatic lead/lag sequencing and automatic alternation of all pumps in order to equalize the amount of usage among the available vacuum pumps. The Premium NFPA Control panel shall include a gateway server card and all features listed below:

- PLC controller and a color touch screen panel which displays the operating status of the unit.
- Building automation communication gateway with BacNet[®] protocol and Web server features. Web servers features include email notifications in case the system is in alarm or has achieved one its maintenance intervals and requires service.
- Ethernet port for connection to BacNet[®] server or direct connection to facility Ethernet for viewing of system operations and status via device IP-address.
- UL508A listed control panel in a NEMA 12 enclosure. The panel door will include: the HMI touch screen, an audible and visual alarms with an acknowledge button, and an HOA switch for each pump.


- Magnetic starters.
- Vacuum transducer for process control.
- Single point power connection.
- Redundant 120Vac control transformers with fused primary and secondary protection.
- System overload trip, high temperature conditions or maintenance intervals for the pump will result in visual and audible alarms.


Available Options

• Internal tank lining for corrosion resistance.



Dimensions										
Model	Dim. A	Dim. B	Dim. C	Inlet	Outlet					
VVOTD0303	36″	53″	90″	1 ½"	1 ½"					
VVOTD0404	36″	54"	89"	1 ½"	1 ½"					
VVOTD0504	37″	52″	90″	1 ½"	1 ½"					





	Oil-less Medical Vacuum Package with Premium Controls												
	Model HP		SCFM @ 19" Hg	NFPA System Capacity ⁽¹⁾	Tank Size (gal)	BTU/Hr ⁽²⁾	dB(A) Level ⁽³⁾	System F.L.A.			System		
		HP						208V	230V	460V	Weight (lbs)		
	VVOTD0303	3	13.5	13.5	80V	7,638	76	17.8	16.4	8.2	830		
	VVOTD0404	5	17	17	120V	12,730	77	24.8	24.8	12.4	1,150		
	VVOTD0504	5	20	20	120V	12,730	76	31.4	29.6	14.8	1,200		

Notes:

- 1 System Capacity is shown with one or more pumps in reserve per NFPA 99.
- 2 BTU/Hr levels are shown with reserve compressor(s) on standby.
- 3 dB(A) levels are shown with one compressor in reserve.
- 4 3 Year Limited Warranty.