




Dimensions					
Model	Dim. A	Dim. B	Dim. C	Inlet	Outlet
MPQ0758	96"	90"	96"	2"	1"
MPQ1008	96"	90"	96"	2"	1 ½"
MPQ1508	96"	90"	96"	2"	1 ½"



Medical Package Quadplex with Premium Controls										
Model	HP	ICFM @ 50 PSIG ^(1,6)	SCFM @ 50 PSIG ^(1,3)	Tank Size (gal)	BTU/Hr ⁽²⁾	dB(A) Level ⁽³⁾	System F.L.A.			System Weight (lbs)
							208V	230V	460V	
MPQ0758	7.5 (4)	87	81.9	240	57,264	81	82	76	39	4,650
MPQ1008	10 (4)	114	108.6	240	76,350	81	111.6	101.2	51.6	4,750
MPQ1508	15 (4)	166.2	152.7	240	114,525	84	166.4	150.8	76.4	4,500

Notes:

- 1 – ICFM and SCFM values @ 50 PSIG are shown with one compressor in reserve per NFPA 99.
- *Powerex recommends using performance ratings in SCFM (Standard Cubic Feet per Minute) when sizing medical air systems.
- *ICFM is standardly defined within the Medical Air Compressor Industry.
- 2 – BTU/Hr levels are shown with reserve pump(s) on standby.
- 3 – dB(A) levels are shown with one pump in reserve per NFPA 99.
- 4 – 3 Year Limited Warranty.