



Claw Vacuum Package Specification

Model	HP	SCFM @ 19" Hg ⁽²⁾	Tank (gal)	BTU/Hr ⁽¹⁾	dB(A) Level ⁽¹⁾	System F.L.A.			Weight (lbs)	Dimensions (in)				
						208V	230V	460V		A	B	C	Inlet	Outlet
CVPT1505-EX4	15 (4)	258 (387)	200	97,539	88	164	142	72	3992 (4788)	99.5	66	83	4	2

- Notes:
- 1 – System Capacity, BTU/Hr, & dB(A) values are shown with one or more compressors in reserve per NFPA 99.
 - 2 – Powerex recommends using performance ratings in SCFM (Standard Cubic Feet per Minute) when sizing medical air systems
 - 3 – Minimum recommended clearance 24" around system; follow OSHA/NEC requirements at control panel.
 - 4 – 3 Year Limited Warranty.