







MSQ15C6xP5

15HP MED OPEN SCROLL QUADPLEX 240GAL DUPLEX DESICCANT

REVISION NOTES	
DESCRIPTION	ECN
REVISED CONTROL PANEL FACE	PXEC0318




- NOTES:
1. APPROXIMATE WEIGHT: 4900 LBS
 2. A MINIMUM OF 24" CLEARANCE IS RECOMMENDED AROUND THE ENTIRE UNIT. CLEARANCE FOR CONTROL PANEL SHALL BE IN ACCORDANCE WITH OSHA/NEC REQUIREMENTS.
 3. SYSTEM DESIGNED IN COMPLIANCE WITH U.S. NFPA '99, LATEST EDITION.
 4. INLET AND OUTLET FLEX HOSES ADD 12" TO CONNECTIONS
1. PESO APROXIMADO: 2223 kg
2. RECOMENDAMOS QUE EL ESPACIO LIBRE ALREDEDOR DEL SISTEMA DEBE SER COMO MÍNIMO 24" (600 mm). EL ESPACIO LIBRE PARA EL TABLERO DE CONTROL DEBERÁ ESTAR DE ACUERDO CON LOS REQUERIMIENTOS DE SEGURIDAD DE LA OSHA/NEC.
3. EL SISTEMA ESTÁ DISEÑADO DE ACUERDO CON LA ÚLTIMA EDICIÓN DE LA NFPA 99.
4. LAS MANGUERAS FLEXIBLES DE ENTRADA Y SALIDA AÑADEN 12" A LAS CONEXIONES

SYMBOL LEGEND		 DRAWINGS TO CONFORM TO CURRENT REVISION OF THE FOLLOWING AND STANDARDS: Y14.5 MULTI AND SECTIONAL, VIEW DWGS Y14.5 DIMENSIONING AND TOLERANCING				
TOLERANCES, UNLESS OTHERWISE SPECIFIED		DO NOT SCALE DRAWING FOR DIMENSIONS				
0.XX ± .05		SOLID EDGE DRAWING		SHEET 1 OF 2	SCALE	1:1
0.XXX ± .005						
± .1		INCH [mm]				
CONFIDENTIAL DISCLOSURE		AUTHOR		REVISION DATE		APPROVALS ON FILE IN PLM
This drawing is the property of POWEREX INC. of MARMON HOLDINGS, INC. and is subject to return on demand. Its contents are confidential and must not be copied or submitted to outside parties for use or examination.		Kumar,Sreehari		07/24/25		
ITEM ID: 761760		 150 PRODUCTION DRIVE, HARRISON OH 45030 (888) 769-7979				
MATERIAL:						

DESCRIPTION						
15HP MED OPEN SCROLL QUADPLEX 240GAL DUPLEX DESICCANT		DRAWING NUMBER		REVISION		SIZE
		MSQ15C6xP5		B		D

15HP MED OPEN SCROLL QUADPLEX 240GAL DUPLEX DESICCANT

POWERX[®]



SYMBOL LEGEND 	DRAWINGS TO CONFORM TO CURRENT REVISION OF THE FOLLOWING ANSI STANDARDS: Y14.3 MULTI AND SECTIONAL VIEW DWGS Y14.5 DIMENSIONING AND TOLERANCING		
TOLERANCES UNLESS OTHERWISE SPECIFIED $0.XX \pm .05$ $0.XXX \pm .005$ $\pm .1$	SOLID EDGE DRAWING	DO NOT SCALE DRAWING FOR DIMENSIONS	
INCH [mm]	SHEET 2 OF 2	SCALE	1:1
* CONFIDENTIAL DISCLOSURE This drawing is the property of POWEREX INC. of WARREN HOLDINGS, INC. and is subject to return on demand. Its contents are confidential and must not be copied or submitted to outside parties for use or examination.	AUTHOR Kumar, Sreehari	REVISION DATE 07/24/25	APPROVALS ON FILE IN PLM
ITEM ID: 761760			
MATERIAL: --	150 PRODUCTION DRIVE, HARRISON OH 45030 (888) 769-7979		
DESCRIPTION 15HP MED OPEN SCROLL QUADPLEX 240GAL DUPLEX	DRAWING NUMBER MSQ15C6xP5 REVISION B SIZE D		