




MSD10C5xP5

10HP MED OPEN SCROLL DUPLEX 200GAL DUPLEX DESICCANT

| REVISION NOTES     |          |
|--------------------|----------|
| DESCRIPTION        | ECN      |
| PRODUCTION RELEASE | PXEC0318 |




- NOTES:
1. APPROXIMATE WEIGHT: 2183 LBS
  2. A MINIMUM OF 24" CLEARANCE IS RECOMMENDED AROUND THE ENTIRE UNIT. CLEARANCE FOR CONTROL PANEL SHALL BE IN ACCORDANCE WITH OSHA/NEC REQUIREMENTS.
  3. SYSTEM DESIGNED IN COMPLIANCE WITH U.S. NFPA '99, LATEST EDITION.
  4. INLET AND OUTLET FLEX HOSES ADD 12" TO CONNECTIONS
1. PESO APROXIMADO: 990 kg
2. RECOMENDAMOS QUE EL ESPACIO LIBRE ALREDEDOR DEL SISTEMA DEBE SER COMO MÍNIMO 24" (600 mm). EL ESPACIO LIBRE PARA EL TABLERO DE CONTROL DEBERÁ ESTAR DE ACUERDO CON LOS REQUERIMIENTOS DE SEGURIDAD DE LA OSHA/NEC.
3. EL SISTEMA ESTÁ DISEÑADO DE ACUERDO CON LA ÚLTIMA EDICIÓN DE LA NFPA 99.
4. LAS MANGUERAS FLEXIBLES DE ENTRADA Y SALIDA AÑADEN 12" A LAS CONEXIONES

| SYMBOL LEGEND                                                                                                                                                                                                                                       |  | DRAWINGS TO CONFORM TO CURRENT REVISION OF THE FOLLOWING AND STANDARDS: |                                     |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--|-------------------------------------------------------------------------|-------------------------------------|
|                                                                                                                                                                                                                                                     |  | Y14.3 MULTI AND SECTIONAL VIEW DWGS                                     |                                     |
|                                                                                                                                                                                                                                                     |  | Y14.5 DIMENSIONING AND TOLERANCING                                      |                                     |
| TOLERANCES, UNLESS OTHERWISE SPECIFIED                                                                                                                                                                                                              |  | SOLID EDGE DRAWING                                                      | DO NOT SCALE DRAWING FOR DIMENSIONS |
| 0.XX ± .05                                                                                                                                                                                                                                          |  |                                                                         |                                     |
| 0.XXX ± .005                                                                                                                                                                                                                                        |  | INCH [mm]                                                               |                                     |
| ± .1                                                                                                                                                                                                                                                |  |                                                                         |                                     |
| CONFIDENTIAL DISCLOSURE<br>This drawing is the property of POWEREX INC. of HARRISON HOLDINGS, INC. and is subject to return on demand. Its contents are confidential and must not be copied or submitted to outside parties for use or examination. |  | AUTHOR                                                                  | REVISION DATE                       |
| DRAWING NUMBER                                                                                                                                                                                                                                      |  | REVISION                                                                | SIZE                                |
| 10HP MED OPEN SCROLL DUPLEX 200GAL DUPLEX DESICCANT                                                                                                                                                                                                 |  | MSD10C5xP5                                                              | B D                                 |



## 10HP MED OPEN SCROLL DUPLEX 200GAL DUPLEX DESICCANT

POWERX<sup>®</sup>



|                                                                                                                                                                                                                                                                   |                       |                                                                                                                                                            |               |
|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|
| SYMBOL LEGEND                                                                                                                                                                                                                                                     |                       | DRAWINGS TO CONFORM TO CURRENT REVISION OF THE FOLLOWING ANSI STANDARDS:<br>Y14.3 MULTI-VIEW AND SECTIONAL VIEW DWGS<br>Y14.5 DIMENSIONING AND TOLERANCING |               |
| TOLERANCES,<br>UNLESS OTHERWISE<br>SPECIFIED<br>0.XX ± .05<br>0.XXX ± .005                                                                                                                                                                                        | SOLID EDGE<br>DRAWING | DO NOT SCALE DRAWING FOR<br>DIMENSIONS                                                                                                                     |               |
| ± .1                                                                                                                                                                                                                                                              | INCH [mm]             | SHEET 2 OF 2                                                                                                                                               | SCALE 1:1     |
| CONFIDENTIAL DISCLOSURE<br>This drawing is the property of POWERS INC.<br>of HARRISON HIGHLANDS.<br>It is subject to return on demand.<br>Its contents are confidential and must be kept<br>and copied or submitted to<br>outside parties for use or examination. |                       | AUTHOR                                                                                                                                                     | REVISION DATE |
|                                                                                                                                                                                                                                                                   |                       | Kumar,Sreehari                                                                                                                                             | 07/10/25      |
| APPROVALS<br>FILE #                                                                                                                                                                                                                                               |                       |                                                                                                                                                            |               |
| ITEM ID:<br>761382                                                                                                                                                                                                                                                |                       | POWEX                                                                                                                                                      |               |
| MATERIAL:<br>--                                                                                                                                                                                                                                                   |                       | 150 PRODUCTION DRIVE, HARRISON OH 45030<br>(888) 769-7079                                                                                                  |               |
| DESCRIPTION<br>10HP MED OPEN SCROLL<br>DUPLEX 200GAL DUPLEX<br>DESICCANT                                                                                                                                                                                          |                       | DRAWING NUMBER<br>MSD10C5Px5                                                                                                                               | REVISION<br>B |
|                                                                                                                                                                                                                                                                   |                       |                                                                                                                                                            | SIZE<br>D     |